

SAFETY DATA SHEET

Revision Date 01-Mar-2017

Revision Number 1

1. IDENTIFICATION OF THE SUBSTANCE/PREPARATION AND OF THE COMPANY/UNDERTAKING

Product identifier			
Product Name	VariKleer Powder		
Product Code(s)	20-3591, 20-3591-002, 20-3591-010		
(M)SDS Number	1344870_A		
Other means of identification			
Synonyms	None		
Recommended use of the chemical	and restrictions on use		
Recommended Use	Laboratory Use Only		
Uses advised against	No information available		
Details of the supplier of the safety data sheet			
Manufacturer	Buehler		
Manufacturer Address	41 Waukegan Rd Lake Bluff, IL 60044 www.buehler.com		
Phone number	+1 847 295 6500		
E-mail Address	custserv@buehler.com		
Emergency telephone number			
Global Access Code: 334545 Americas: +1 760 476 3962 Middle East/Africa: +1 760 476 3959	Asia Pacific: +1 760 476 3960 Europe: +1 760 476 3961		

2. HAZARDS IDENTIFICATION

Classification

This chemical is considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200).

Skin sensitization

Category 1

GHS Label elements, including precautionary statements

Emergency Overview



Signal word	Warning	
Hazard Statements May cause an allergic skin reaction		
Appearance Colorless	Physical state Powder(s)	Odor Characteristic

Precautionary Statements - Prevention

Avoid breathing dust/fume/gas/mist/vapors/spray Contaminated work clothing should not be allowed out of the workplace Wear protective gloves

Precautionary Statements - Response None

Skin

IF ON SKIN: Wash with plenty of soap and water If skin irritation or rash occurs: Get medical advice/attention

Precautionary Statements - Storage

None

Precautionary Statements - Disposal

Dispose of contents/container to an approved waste disposal plant

Hazards not otherwise classified (HNOC)

Not applicable

Other information

No information available

Interactions with Other Chemicals

No information available.

3. COMPOSITION/INFORMATION ON INGREDIENTS

Chemical name	CAS No	Weight-%	Trade Secret
Polymethyl methacrylate	9011-14-7	90 - 100%	*
Methyl methacrylate	80-62-6	0 - 1%	*
Benzoyl peroxide	94-36-0	0 - 1%	*

*The exact percentage (concentration) of composition has been withheld as a trade secret



4. FIRST AID MEASURES

First aid measures

General Advice	Show this safety data sheet to the doctor in attendance.	
Eye contact	Rinse thoroughly with plenty of water, also under the eyelids. If symptoms persist, call a physician.	
Skin contact	Wash with soap and water. May cause an allergic skin reaction. In the case of skin irritation or allergic reactions see a physician.	
Inhalation	Remove to fresh air.	
Ingestion	Rinse mouth immediately and drink plenty of water. Never give anything by mouth to an unconscious person.	
Most important symptoms and effe	cts, both acute and delayed	
Most Important Symptoms and Effects	Itching. Rashes. Hives.	
Indication of any immediate medica	dical attention and special treatment needed	
Notes to Physician	May cause sensitization in susceptible persons. Treat symptomatically.	

5. FIRE-FIGHTING MEASURES

Suitable Extinguishing Media

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

Unsuitable extinguishing media

CAUTION: Use of water spray when fighting fire may be inefficient.

Specific hazards arising from the chemical

Product is or contains a sensitizer. May cause sensitization by skin contact.

Uniform Fire Code Sensitizer: Solid

Explosion Data Sensitivity to Mechanical Impact None.

Sensitivity to Static Discharge None.

Protective equipment and precautions for firefighters

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.



6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emergency procedures

Personal precautions	Avoid contact with skin, eyes or clothing. Use personal protective equipment as required.
Environmental precautions	
Environmental precautions	See Section 12 for additional Ecological Information.
Methods and material for containm	ent and cleaning up
Methods for containment	Prevent further leakage or spillage if safe to do so.
Methods for cleaning up	Pick up and transfer to properly labeled containers.

7. HANDLING AND STORAGE

Precautions for safe handling

Handling Handle in accordance with good industrial hygiene and safety practice. Avoid contact with skin, eyes or clothing. Do not eat, drink or smoke when using this product. Take off contaminated clothing and wash before reuse.

Conditions for safe storage, including any incompatibilities

Storage	Keep containers tightly closed in a dry, cool and well-ventilated place. Store locked up. Keep out of the reach of children.
Incompatible Products	None known based on information supplied.

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Control parameters

Exposure Guidelines

This product, as supplied, does not contain any hazardous materials with occupational exposure limits established by the region specific regulatory bodies

Chemical name	ACGIH TLV	OSHA PEL	NIOSH IDLH
Methyl methacrylate	STEL: 100 ppm	TWA: 100 ppm	IDLH: 1000 ppm
80-62-6	TWA: 50 ppm	TWA: 410 mg/m ³	TWA: 100 ppm
		(vacated) TWA: 100 ppm	TWA: 410 mg/m ³
		(vacated) TWA: 410 mg/m ³	_
Benzoyl peroxide	TWA: 5 mg/m ³	TWA: 5 mg/m ³	IDLH: 1500 mg/m ³
94-36-0		(vacated) TWA: 5 mg/m ³	TWA: 5 mg/m ³

ACGIH TLV: American Conference of Governmental Industrial Hygienists - Threshold Limit Value OSHA PEL: Occupational Safety and Health Administration - Permissible Exposure Limits NIOSH IDLH Immediately Dangerous to Life or Health

Other Exposure Guidelines

Vacated limits revoked by the Court of Appeals decision in AFL-CIO v. OSHA, 965 F.2d 962 (11th Cir., 1992)

Appropriate engineering controls



Engineering Measures	Showers Eyewash stations Ventilation systems	
Individual protection measures, such as personal protective equipment		
Eye/face protection	Wear safety glasses with side shields (or goggles).	
Skin and body protection	Wear protective gloves and protective clothing.	
Respiratory protection	If exposure limits are exceeded or irritation is experienced, NIOSH/MSHA approved respiratory protection should be worn. Positive-pressure supplied air respirators may be required for high airborne contaminant concentrations. Respiratory protection must be provided in accordance with current local regulations.	
Hygiene Measures	Handle in accordance with good industrial hygiene and safety practice.	
9. PHYSICAL AND CHEMICAL PROPERTIES		

Physical and Chemical Properties

Physical state	Powder(s)		
Appearance	Colorless	Odor	Characteristic
Color	No information available	Odor Threshold	No information available
Property_	Values	Remarks Method	
pH	No data available	None known	
Melting / freezing point	No data available	None known	
Boiling point / boiling range	No data available	None known	
Flash Point	>250°C	None known	
Evaporation Rate	No data available	None known	
Flammability (solid, gas)	No data available	None known	
Flammability Limit in Air			
Upper flammability limit	No data available		
Lower flammability limit	No data available		
Vapor pressure	No data available	None known	
Vapor density	No data available	None known	
Specific Gravity	No data available	None known	
Water Solubility	Insoluble		
Solubility in other solvents	No data available	None known	
Partition coefficient: n-octanol/wat	erNo data available	None known	
Autoignition temperature	No data available	None known	
Decomposition temperature	No data available	None known	
Kinematic viscosity	No data available	None known	
Dynamic viscosity	No data available	None known	
Explosive properties	No data available		
Oxidizing properties	No data available		
Other Information			
Softening Point	No data available		
VOC Content (%)	No data available		
Particle Size	No data available		
Particle Size Distribution			



10. STABILITY AND REACTIVITY

Reactivity

No data available.

Chemical stability

Stable under recommended storage conditions.

Possibility of Hazardous Reactions

None under normal processing.

Conditions to avoid

None known based on information supplied.

Incompatible materials

None known based on information supplied.

Hazardous Decomposition Products

None known based on information supplied.

11. TOXICOLOGICAL INFORMATION

Information on likely routes of exposure

Product Information

Inhalation	Specific test data for the substance or mixture is not available.
Eye contact	Specific test data for the substance or mixture is not available.
Skin contact	Specific test data for the substance or mixture is not available.
Ingestion	Specific test data for the substance or mixture is not available.

Component Information

Chemical name	Oral LD50	Dermal LD50	Inhalation LC50
Methyl methacrylate 80-62-6	= 7900 mg/kg(Rat)= 7872 mg/kg(Rat)	> 5 g/kg(Rabbit)	= 4632 ppm(Rat)4 h
Benzoyl peroxide 94-36-0	= 7710 mg/kg(Rat)= 6400 mg/kg(Rat)	-	-

Information on toxicological effects

Symptoms Itching. Rashes. Hives.

Delayed and immediate effects as well as chronic effects from short and long-term exposure

Sensitization	May cause sensitization in susceptible persons. May cause sensitization by skin contact.
Mutagenic Effects	No information available.
Carcinogenicity	Contains no ingredient listed as a carcinogen.



Chemical name	ACGIH	IARC	NTP	OSHA
Polymethyl methacrylate 9011-14-7		Group 3		
Methyl methacrylate		Group 3		
80-62-6 Benzoyl peroxide		Group 3		
94-36-0				

IARC (International Agency for Research on Cancer) Group 3 - Not Classifiable as to Carcinogenicity in Humans

Reproductive toxicity No information available.	
STOT - single exposure No information available.	
STOT - repeated exposure No information available.	
Chronic Toxicity	Carcinogenic potential is unknown.
Target Organ Effects	Skin.
Aspiration Hazard	No information available.

Numerical measures of toxicity Product Information

The following values are calculated based on chapter 3.1 of the GHS document Not applicable



12. ECOLOGICAL INFORMATION

Ecotoxicity

The environmental impact of this product has not been fully investigated.

Chemical name	Toxicity to Algae	Toxicity to Fish	Toxicity to Microorganisms	Daphnia Magna (Water Flea)
Methyl methacrylate 80-62-6	96h EC50: = 170 mg/L (Pseudokirchneriella subcapitata)	96h LC50: 243 - 275 mg/L (Pimephales promelas) 96h LC50: 125.5 - 190.7 mg/L (Pimephales promelas) 96h LC50: 153.9 - 341.8 mg/L (Lepomis macrochirus) 96h LC50: > 79 mg/L (Oncorhynchus mykiss) 96h LC50: 170 - 206 mg/L (Lepomis macrochirus) 96h LC50: 326.4 - 426.9 mg/L (Poecilia reticulata)		48h EC50: = 69 mg/L

Persistence and Degradability

No information available.

Bioaccumulation

Chemical name	Log Pow
Methyl methacrylate	0.7
80-62-6	

Other adverse effects

No information available.

13. DISPOSAL CONSIDERATIONS

Waste treatment methods

Disposal methods	This material, as supplied, is a hazardous waste according to federal regulations (40 CFR 261).
Contaminated Packaging	Dispose of contents/containers in accordance with local regulations.
US EPA Waste Number	U162

This product contains one or more substances that are listed with the State of California as a hazardous waste.

Chemical name	California Hazardous Waste
Methyl methacrylate	Toxic
80-62-6	Ignitable
Benzoyl peroxide	Toxic
94-36-0	Ignitable
	Reactive



14. TRANSPORT INFORMATION

ADN	NOT REGULATED
ADR	NOT REGULATED
<u>RID</u>	NOT REGULATED
IMDG/IMO	NOT REGULATED
IATA_ Proper Shipping Name	NOT REGULATED NON REGULATED
ICAO	NOT REGULATED
MEX	NOT REGULATED
TDG	NOT REGULATED
<u>DOT</u> Proper Shipping Name Hazard Class	NOT REGULATED NON REGULATED N/A

15. REGULATORY INFORMATION

International Inventories

TSCA	Complies
DSL	All components are listed either on the DSL or NDSL.

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List

US Federal Regulations

SARA 313

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product contains a chemical or chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372

Chemical name	CAS No	Weight-%	SARA 313 - Threshold Values %
Methyl methacrylate - 80-62-6	80-62-6	0 - 1%	1.0
Benzoyl peroxide - 94-36-0	94-36-0	0 - 1%	1.0
SARA 311/312 Hazard Categories			
Acute Health Hazard	Yes		
Chronic Health Hazard	No		
Fire Hazard	No		
Sudden release of pressure hazard	No		
Reactive Hazard	No		

<u>CWA (Clean Water Act)</u> This product contains the following substances which are regulated pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42)

Chemical name	CWA - Reportable Quantities	CWA - Toxic Pollutants	CWA - Priority Pollutants	CWA - Hazardous Substances
Methyl methacrylate 80-62-6	1000 lb			Х
CERCLA				



This material, as supplied, contains one or more substances regulated as a hazardous substance under the Comprehensive Environmental Response Compensation and Liability Act (CERCLA) (40 CFR 302)

Chemical name	Hazardous Substances RQs	Extremely Hazardous Substances RQs	RQ
Methyl methacrylate 80-62-6	1000 lb		RQ 1000 lb final RQ RQ 454 kg final RQ

US State Regulations

California Proposition 65

This product does not contain any Proposition 65 chemicals.

U.S. State Right-to-Know Regulations

Chemical name	New Jersey	Massachusetts	Pennsylvania	Rhode Island	Illinois
Methyl methacrylate 80-62-6	Х	Х	Х	Х	Х
Benzoyl peroxide 94-36-0	Х	Х	Х	Х	

International Regulations

Component	Carcinogen Status	Exposure Limits
Methyl methacrylate		Mexico: TWA 100 ppm
80-62-6 (0 - 1%)		Mexico: TWA 410 mg/m ³
		Mexico: STEL 125 ppm
		Mexico: STEL 510 mg/m ³
Benzoyl peroxide		Mexico: TWA 5 mg/m ³
94-36-0 (0 - 1%)		

Canada

WHMIS Hazard Class

Not determined

16. OTHER INFORMATION					
NFPA	Health Hazards 2	Flammability 0	Instability 0	Physical and Chemical Hazards -	
HMIS	Health Hazards 2	Flammability 0	Physical Hazard 0	Personal Protection X	
Prepared By	Product Stewardship 23 British American Blvd. Latham, NY 12110 1-800-572-6501				
Issuing Date	01-Mar-2017				
Revision Date	01-Mar-2017				
Revision Note	No information available				

Disclaimer

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text





The supplier identified below generated this SDS using the UL SDS template. UL did not test, certify, or approve the substance described in this SDS, and all information in this SDS was provided by the supplier or was reproduced from publically available regulatory data sources. UL makes no representations or warranties regarding the completeness or accuracy of the information in this SDS and disclaims all liability in connection with the use of this information or the substance described in this SDS. The layout, appearance and format of this SDS is © 2014 UL LLC. All rights reserved.

End of Safety Data Sheet

