Safety Data Sheet



Section 1: Identification of the Substance/Mixture and of the Company/Undertaking

1.1 Product identifier

Product Name Levigated Alumina / SDS# 9104895

Product Code • 40-6435-016, 40-6435-080

1.2 Relevant identified uses of the substance or mixture and uses advised against

Relevant identified use(s) • Polishing Compound

1.3 Details of the supplier of the safety data sheet

Manufacturer
 BUEHLER, a division of Illinios Tool Works Inc.

41 Waukegan Road Lake Bluff, IL 60044 United States

https://www.buehler.com

Telephone (Technical) • 847-295-6500

1.4 Emergency telephone number

Manufacturer • 800-424-9300 - CHEMTREC

Section 2: Hazards Identification

EU/EEC

According to: Regulation (EC) No 1272/2008 (CLP)/REACH 1907/2006 [amended by 453/2010]

According to: EU Directive 67/548/EEC (DSD) or 1999/45/EC (DPD)

2.1 Classification of the substance or mixture

CLP • Specific Target Organ Toxicity Repeated Exposure 2 - H373

DSD/DPD • Harmful (Xn)

R48/20

2.2 Label Elements

CLP

WARNING



Hazard statements • H373 - May cause damage to organs through prolonged or repeated exposure.

Precautionary statements

Prevention • P260 - Do not breathe dust.

Response • P314 - Get medical advice/attention if you feel unwell.

Storage/Disposal • P501 - Dispose of content and/or container in accordance with local, regional,

national, and/or international regulations.

DSD/DPD



Risk phrases • R48/20 - Harmful: danger of serious damage to health by prolonged exposure through

inhalation.

2.3 Other Hazards

• According to Regulation (EC) No. 1272/2008 (CLP) this material is considered

hazardous.

• According to European Directive 1999/45/EC this material is considered dangerous.

UN GHS

According to: UN Globally Harmonized System of Classification and Labelling of Chemicals (GHS)

2.1 Classification of the substance or mixture

UN GHS

Specific Target Organ Toxicity Repeated Exposure 2

2.2 Label elements

UN GHS

WARNING



Hazard statements • May cause damage to organs through prolonged or repeated exposure.

Precautionary statements

Prevention • Do not breathe dust.

Response • Get medical advice/attention if you feel unwell.

Storage/Disposal • Dispose of content and/or container in accordance with local, regional, national, and/or

international regulations.

2.3 Other hazards

According to the Globally Harmonized System for Classification and Labeling (GHS)

this product is considered hazardous.

United States (US)

According to: OSHA 29 CFR 1910.1200 HCS

2.1 Classification of the substance or mixture

OSHA HCS 2012

Specific Target Organ Toxicity Repeated Exposure 2

2.2 Label elements

OSHA HCS 2012

WARNING



Hazard statements • May cause damage to organs through prolonged or repeated exposure.

Precautionary statements

Prevention • Do not breathe dust.

Response • Get medical advice/attention if you feel unwell.

Storage/Disposal • Dispose of content and/or container in accordance with local, regional, national, and/or

international regulations.

2.3 Other hazards

• Under United States Regulations (29 CFR 1910.1200 - Hazard Communication

Standard), this product is considered hazardous.

Canada

According to: WHMIS

2.1 Classification of the substance or mixture

WHMIS
 Other Toxic Effects - D2B

2.2 Label elements

WHMIS



· Other Toxic Effects - D2B

2.3 Other hazards

WHMIS

• In Canada, the product mentioned above is considered hazardous under the Workplace Hazardous Materials Information System (WHMIS).

Section 3 - Composition/Information on Ingredients

3.1 Substances

· Material does not meet the criteria of a substance.

3.2 Mixtures

Composition							
Chemical Identifiers		%	LD50/LC50	Classifications According to Regulation/Directive	Comments		
Aluminum oxide	CAS:1344-28-1 EC Number:215- 691-6	70% TO 100%	Inhalation-Rat LC50 • 0.2 mg/L 5 Hour(s) 28 Week(s)	UN GHS: STOT RE 2 (Lungs, Inhl); EU DSD/DPD: Xn, R48/20; EU CLP: STOT RE 2, H373 (Lungs, Inhl); OSHA HCS 2012: STOT RE 2 (Lungs, Inhl);	NDA		

Section 4 - First Aid Measures

4.1 Description of first aid measures

Inhalation

 Move victim to fresh air. Give artificial respiration if victim is not breathing. Administer oxygen if breathing is difficult. If signs/symptoms continue, get medical attention.

Skin

In case of contact with substance, immediately flush skin with running water for at least 20 minutes. If irritation develops and persists, get medical attention.

Eye

In case of contact with substance, immediately flush eyes with running water for at least 20 minutes. If eye irritation persists: Get medical advice/attention.

Ingestion

Do NOT induce vomiting. Give victim a glass of water or milk. Never give anything by mouth to an unconscious person. Obtain medical attention immediately if ingested.

4.2 Most important symptoms and effects, both acute and delayed

· Refer to Section 11 - Toxicological Information.

4.3 Indication of any immediate medical attention and special treatment needed

Notes to Physician

 All treatments should be based on observed signs and symptoms of distress in the patient. Consideration should be given to the possibility that overexposure to materials other than this product may have occurred.

Section 5 - Firefighting Measures

5.1 Extinguishing media

Suitable Extinguishing Media • In case of fire use media as appropriate for surrounding fire.

Unsuitable Extinguishing Media

No data available.

5.2 Special hazards arising from the substance or mixture

Unusual Fire and Explosion

Hazards

· None known.

Hazardous Combustion

Products

No data available.

5.3 Advice for firefighters

Wear positive pressure self-contained breathing apparatus (SCBA). Structural firefighters' protective clothing will only provide limited protection.

Section 6 - Accidental Release Measures

6.1 Personal precautions, protective equipment and emergency procedures

Personal Precautions

 Ventilate enclosed areas. Do not walk through spilled material. Wear appropriate personal protective equipment, avoid direct contact.

Emergency Procedures

As an immediate precautionary measure, isolate spill or leak area for at least 25 meters (75 feet) in all directions. If tank, rail car or tank truck is involved in a fire, ISOLATE for 800 meters (1/2 mile) in all directions; also, consider initial evacuation for 800 meters (1/2 mile) in all directions. Keep unauthorized personnel away.

6.2 Environmental precautions

Avoid run off to waterways and sewers.

6.3 Methods and material for containment and cleaning up

Containment/Clean-up Measures

Avoid generating dust.

SMALL DRY SPILLS: With clean shovel place material into clean, dry container and cover loosely; move containers from spill area.

LARGE SPILLS: Cover powder spill with plastic sheet or tarp to minimize spreading.

6.4 Reference to other sections

 Refer to Section 8 - Exposure Controls/Personal Protection and Section 13 - Disposal Considerations.

Section 7 - Handling and Storage

7.1 Precautions for safe handling

Handling

 Use only with adequate ventilation. Wear appropriate personal protective equipment, avoid direct contact. Do not breathe dust. Avoid contact with skin, eyes, and clothing. Wash thoroughly with soap and water after handling and before eating, drinking, or using tobacco.

7.2 Conditions for safe storage, including any incompatibilities

Storage

 Keep container tightly closed. Store in a cool, dry, well-ventilated place. Do not store above 24°C (75.2°F).

7.3 Specific end use(s)

Refer to Section 1.2 - Relevant identified uses.

Section 8 - Exposure Controls/Personal Protection

8.1 Control parameters

Exposure Limits/Guidelines			
	Result	ACGIH	OSHA
Aluminum oxide	TWAs		15 mg/m3 TWA (total dust); 5 mg/m3 TWA (respirable
(1344-28-1)		as Aluminum insoluble compounds	fraction)

8.2 Exposure controls

Engineering Measures/Controls

Good general ventilation should be used. Ventilation rates should be matched to
conditions. If applicable, use process enclosures, local exhaust ventilation, or other
engineering controls to maintain airborne levels below recommended exposure limits.
If exposure limits have not been established, maintain airborne levels to an acceptable
level. Ensure that dust handling systems (such as exhaust ducts, dust collectors,
vessels and processing equipment) are designed in a manner to prevent the escape of
dust into the work area (i.e., there is not leakage from the equipment).

Personal Protective Equipment

Respiratory

For limited exposure use an N95 dust mask. For prolonged exposure use an airpurifying respirator with high efficiency particulate air (HEPA) filters. Follow the OSHA
respirator regulations found in 29 CFR 1910.134 or European Standard EN 149. Use a
NIOSH/MSHA or European Standard EN 149 approved respirator if exposure limits are
exceeded or symptoms are experienced.

Eye/Face

Wear safety goggles.

Skin/Body

Wear appropriate gloves. Wear long sleeves and/or protective coveralls.

Environmental Exposure Controls

Controls should be engineered to prevent release to the environment, including
procedures to prevent spills, atmospheric release and release to waterways. Follow
best practice for site management and disposal of waste.

Key to abbreviations

ACGIH = American Conference of Governmental Industrial Hygiene

OSHA = Occupational Safety and Health Administration

TWA = Time-Weighted Averages are based on 8h/day, 40h/week exposures

Section 9 - Physical and Chemical Properties

9.1 Information on Physical and Chemical Properties

Material Description						
Physical Form	Solid	Appearance/Description	White powder with no odor.			
Color	White	Odor	None			
Odor Threshold	Data lacking					
General Properties						
Boiling Point	Data lacking	Melting Point/Freezing Point	Data lacking			
Decomposition Temperature	Data lacking	рН	Data lacking			
Specific Gravity/Relative Density	= 3.8425 Water=1	Water Solubility	Insoluble			
Viscosity	Data lacking	Explosive Properties	Data lacking			
Oxidizing Properties:	Data lacking					
Volatility						
Vapor Pressure	Data lacking	Vapor Density	< 1 Air=1			
Evaporation Rate	< 1 n-Butyl Acetate = 1					
Flammability						
Flash Point	Data lacking	UEL	Data lacking			
LEL	Data lacking	Autoignition	Data lacking			
Flammability (solid, gas)	Data lacking					
Environmental						
Octanol/Water Partition coefficient	Data lacking					

9.2 Other Information

No additional physical and chemical parameters noted.

Section 10: Stability and Reactivity

10.1 Reactivity

· No dangerous reaction known under conditions of normal use.

10.2 Chemical stability

Stable under normal temperatures and pressures.

10.3 Possibility of hazardous reactions

· Hazardous polymerization not indicated.

10.4 Conditions to avoid

· Avoid generating dust.

10.5 Incompatible materials

 Reactive with oxidizing agents, acids. Chlorine Trifluoride reacts violently with Aluminum Oxide producing a flame. Ethylene oxide may polymerize violently when in contact with highly catalytic surfaces such as pure Aluminum Oxide. Reacts with hot chlorinated rubber.

10.6 Hazardous decomposition products

· No data available.

Section 11 - Toxicological Information

11.1 Information on toxicological effects

	Components					
Loxide (70% TOL	1344- 28-1	Multi-dose Toxicity: Inhalation-Rat TCLo • 200 mg/m³ 5 Hour(s) 28 Week(s)-Intermittent; Lungs, Thorax, or Respiration:Structural or functional change in trachea or bronchi; Lungs, Thorax, or Respiration:Chronic pulmonary edema; Related to Chronic Data:Death in the Other Multiple Dose data type field; Tumorigen / Carcinogen: Implant-Rat • 200 mg/kg; Tumorigenic:Equivocal tumorigenic agent by RTECS criteria; Tumorigenic:Tumors at site of application; Implant-Rat TDLo • 200 mg/kg; Tumorigenic:Neoplastic by RTECS criteria; Tumorigenic:Tumors at site of application; Intrapleural-Rat TDLo • 90 mg/kg; Tumorigenic:Equivocal tumorigenic agent by RTECS criteria; Lungs, Thorax, or Respiration:Tumors				

GHS Properties	Classification
Respiratory sensitization	EU/CLP • Data lacking OSHA HCS 2012 • Data lacking UN GHS • Data lacking
Serious eye damage/Irritation	EU/CLP • Data lacking OSHA HCS 2012 • Data lacking UN GHS • Data lacking
Acute toxicity	EU/CLP • Data lacking OSHA HCS 2012 • Data lacking UN GHS • Data lacking
Aspiration Hazard	EU/CLP • Data lacking OSHA HCS 2012 • Data lacking UN GHS • Data lacking
Carcinogenicity	EU/CLP • Data lacking OSHA HCS 2012 • Data lacking UN GHS • Data lacking
Skin corrosion/Irritation	EU/CLP • Data lacking OSHA HCS 2012 • Data lacking UN GHS • Data lacking
Skin sensitization	EU/CLP • Data lacking OSHA HCS 2012 • Data lacking UN GHS • Data lacking
STOT-RE	EU/CLP • Specific Target Organ Toxicity Repeated Exposure 2 OSHA HCS 2012 • Specific Target Organ Toxicity Repeated Exposure 2 UN GHS • Specific Target Organ Toxicity Repeated Exposure 2
STOT-SE	EU/CLP • Data lacking OSHA HCS 2012 • Data lacking UN GHS • Data lacking
Toxicity for Reproduction	EU/CLP • Data lacking OSHA HCS 2012 • Data lacking UN GHS • Data lacking
Germ Cell Mutagenicity	EU/CLP • Data lacking OSHA HCS 2012 • Data lacking UN GHS • Data lacking

Potential Health Effects Inhalation

Acute (Immediate)

 Exposure to dust may cause irritation. Processes such as cutting, grinding, crushing, or impact may result in generation of excessive amounts of airborne dusts in the workplace. Nuisance dust may affect the lungs but reactions are typically reversible.

Chronic (Delayed)

 Repeated and prolonged exposure to dust may cause lung effects including pneumoconiosis.

Skin

Acute (Immediate)

Exposure to dust may cause mechanical irritation.

Chronic (Delayed)

· No data available.

Eye

Acute (Immediate)

 Exposure to dust may cause mechanical irritation. Excessive concentrations of nuisance dust in the workplace may reduce visibility and may cause unpleasant deposits in eyes.

Chronic (Delayed)

· No data available.

Ingestion

Acute (Immediate)

• Excessive concentrations of nuisance dust in the workplace may cause mechanical irritation to mucous membranes.

Chronic (Delayed)

No data available.

Key to abbreviations

TC = Toxic Concentration

TD = Toxic Dose

Section 12 - Ecological Information

12.1 Toxicity

· Material data lacking.

12.2 Persistence and degradability

 The product itself and its products of degradation are not toxic. Possibly hazardous short term degradation products are not likely. However, long term degradation products may arise.

12.3 Bioaccumulative potential

Material data lacking.

12.4 Mobility in Soil

· Material data lacking.

12.5 Results of PBT and vPvB assessment

No PBT and vPvB assessment has been conducted.

12.6 Other adverse effects

No studies have been found.

Section 13 - Disposal Considerations

13.1 Waste treatment methods

Product waste

• Dispose of content and/or container in accordance with local, regional, national, and/or international regulations.

Packaging waste

 Dispose of content and/or container in accordance with local, regional, national, and/or international regulations.

Section 14 - Transport Information

Preparation Date: 03/April/2012 Revision Date: 28/October/2015

	14.1 UN number	14.2 UN proper shipping name	14.3 Transport hazard class(es)	14.4 Packing group	14.5 Environmental hazards
DOT	NDA	Not Regulated	NDA	NDA	NDA
TDG	NDA	Not Regulated	NDA	NDA	NDA
IATA/ICAO	NDA	Not Regulated	NDA	NDA	NDA

14.6 Special precautions for • None specified.

14.7 Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Data lacking.

Section 15 - Regulatory Information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

SARA Hazard Classifications • Chronic

			Inventory			
Component	CAS	Canada DSL	Canada NDSL	EU EINECS	EU ELNICS	TSCA
Aluminum oxide	1344-28-1	Yes	No	Yes	No	Yes

Canada

-Labor		
Canada - W	MIS - Classifications of Substances	

Aluminum oxide

1344-28-1

Uncontrolled product according to WHMIS classification criteria

Canada - WHMIS - Ingredient Disclosure List

· Aluminum oxide

1344-28-1

Environment

Canada - CEPA - Priority Substances List

Aluminum oxide

1344-28-1 Not Listed

1 %

United States

Labor

U.S. - OSHA - Process Safety Management - Highly Hazardous Chemicals

· Aluminum oxide

Not Listed 1344-28-1

U.S. - OSHA - Specifically Regulated Chemicals

· Aluminum oxide

1344-28-1 Not Listed

Environment

U.S. - CAA (Clean Air Act) - 1990 Hazardous Air Pollutants

· Aluminum oxide

1344-28-1 Not Listed

U.S. - CERCLA/SARA - Hazardous Substances and their Reportable Quantities

· Aluminum oxide

1344-28-1 Not Listed

Preparation Date: 03/April/2012 Revision Date: 28/October/2015

U.S CERCLA/SARA - Radionuclides and Their Reportable Quantities • Aluminum oxide	1344-28-1	Not Listed	
U.S CERCLA/SARA - Section 302 Extremely Hazardous Substances EPCRA RQs • Aluminum oxide	1344-28-1	Not Listed	
U.S CERCLA/SARA - Section 302 Extremely Hazardous Substances TPQs • Aluminum oxide	1344-28-1	Not Listed	
U.S CERCLA/SARA - Section 313 - Emission Reporting		1.0 % de minimis	
Aluminum oxide	1344-28-1	concentration (fibrous forms)	
U.S CERCLA/SARA - Section 313 - PBT Chemical Listing • Aluminum oxide	1344-28-1	Not Listed	

United States - California

- Environment		
U.S California - Proposition 65 - Carcinogens ListAluminum oxide	1344-28-1	Not Listed
U.S California - Proposition 65 - Developmental Toxicity • Aluminum oxide	1344-28-1	Not Listed
U.S California - Proposition 65 - Maximum Allowable Dose Levels (MADL) • Aluminum oxide	1344-28-1	Not Listed
U.S California - Proposition 65 - No Significant Risk Levels (NSRL) • Aluminum oxide	1344-28-1	Not Listed
U.S California - Proposition 65 - Reproductive Toxicity - Female • Aluminum oxide	1344-28-1	Not Listed
U.S California - Proposition 65 - Reproductive Toxicity - Male • Aluminum oxide	1344-28-1	Not Listed

15.2 Chemical Safety Assessment

No Chemical Safety Assessment has been carried out.

Section 16 - Other Information

Revision Date

· 28/October/2015

Preparation Date

• 03/April/2012

Disclaimer/Statement of Liability

• To the best of our knowledge, the information contained in this SDS is accurate or is obtained from sources believed to be accurate. However, no liability, expressed or implied, is assumed for the accuracy or completeness of the information contained herein. Buyer assumes liability in its use of the material.

Key to abbreviations

NDA = No data available